

Create Account			
Test Case ID	Test Scenario	Preconditions	Test steps
TC_001	Create Account with Valid Input	User is on the registration page	1. Enter valid first name (e.g., "John"). 2. Enter valid last name (e.g., "Doe"). 3. Enter valid username (e.g., "johndoe"). 4. Enter valid password (e.g., "Password123"). 5. Confirm password (e.g., "Password123"). 6. Enter valid email (e.g., "johndoe@example.com"). 7. Select valid date of birth from the drop-down calendar (e.g., "01/01/1990"). 8. Click "Create Account".
TC_002	Create Account with Invalid Email	User is on the registration page.	1. Enter valid first name (e.g., "John"). 2. Enter valid last name (e.g., "Doe"). 3. Enter valid username (e.g., "johndoe"). 4. Enter valid password (e.g., "Password123"). 5. Confirm password (e.g., "Password123"). 6. Enter invalid email (e.g., "johndoe"). 7. Select valid date of birth from the drop-down calendar (e.g., "01/01/1990"). 8. Click "Create Account".
TC_003	Create Account with Mismatched Passwords	User is on the registration page.	1. Enter valid first name (e.g., "John"). 2. Enter valid last name (e.g., "Doe"). 3. Enter valid username (e.g., "johndoe"). 4. Enter valid password (e.g., "Password123"). 5. Confirm password (e.g., "Password456"). 6. Enter valid email (e.g., "johndoe@example.com"). 7. Select valid date of birth from the drop-down calendar (e.g., "01/01/1990"). 8. Click "Create Account".
Edit Mealtimes			
TC_001	Edit Mealtimes with Valid Input	User is logged in. Carb-to-unit ratio is stored (e.g., 10).	1. Navigate to "Edit Mealtimes & Algorithm". 2. Select new mealtime from the drop-down menu (e.g., "08:00 AM"). 3. Enter valid carb-to-unit ratio (e.g., "15"). 4. Click "Save".
TC_002	Edit Mealtimes with Missing Carb-to-Unit Ratio	User is logged in. Carb-to-unit ratio is stored (e.g., 10).	1. Navigate to "Edit Mealtimes & Algorithm". 2. Select new mealtime from the drop-down menu (e.g., "08:00 AM"). 3. Leave carb-to-unit ratio field empty. 4. Click "Save".
TC_003	Edit Mealtimes with Default Mealtime (00:00 AM/PM)	User is logged in. Carb-to-unit ratio is stored (e.g., 10).	1. Navigate to "Edit Mealtimes & Algorithm". 2. Leave mealtime as the default value (00:00 AM/PM). 3. Enter valid carb-to-unit ratio (e.g., "15"). 4. Click "Save".
TC_004	Edit Mealtimes with Multiple Mealtimes	User is logged in. Carb-to-unit ratio is stored (e.g., 10).	1. Navigate to "Edit Mealtimes & Algorithm". 2. Select new mealtime for breakfast (e.g., "08:00 AM"). 3. Enter valid carb-to-unit ratio for breakfast (e.g., "15"). 4. Select new mealtime for lunch (e.g., "12:00 PM"). 5. Enter valid carb-to-unit ratio for lunch (e.g., "10"). 6. Click "Save".
TC_005	Edit Mealtimes and Navigate Away Without Saving	User is logged in. Carb-to-unit ratio is stored (e.g., 10)	1. Navigate to "Edit Mealtimes & Algorithm". 2. Select new mealtime from the drop-down menu (e.g., "08:00 AM"). 3. Enter valid carb-to-unit ratio (e.g., "15"). 4. Navigate away from the edit screen without clicking "Save".
Authentication			
TC_Auth_001	System authenticates a user with valid credentials.	User credentials exist in the database and user is not already logged in.	1. Navigate to login screen. 2. Enter valid username and password. 3. Click "Login".
TC_Auth_002	System rejects authentication with a username not in the database.	Invalid username.	1. Enter a non-existent username. 2. Enter any password. 3. Click "Login".
TC_Auth_003	System rejects authentication with incorrect password.	Username exists but password is incorrect.	1. Enter a valid username. 2. Enter wrong password. 3. Click "Login".
TC_Auth_004	System prevents duplicate login for same credentials if already logged in.	User session already exists.	1. Attempt login from a second device/tab.
TC_Auth_005	System should not authenticate if username/password fields are empty.	Login page is loaded.	1. Leave one or both fields empty. 2. Click "Login".
TC_Auth_006	Verify login immediately after successful sign-up.	User completes sign-up.	1. Sign up with new username/password. 2. Attempt to log in with those credentials.
TC_Auth_007	On successful authentication, user adherence data should load.	User has adherence records stored.	1. Login with valid credentials.
Access Data			
TC_001	User successfully accesses data	User is on the main screen User has data to access	1. Navigate to main screen 2. Press access data button 3. Scroll to find relevant data
TC_002	User does not have data to access	User is on the main screen User has not collected any data	1. Navigate to main screen 2. Press access data button 3. Error message is displayed
TC_003	User exits data display	User is accessing data	1. Press back button

Expected result	Postconditions	
Account is created successfully. User is redirected to the login screen.	User credentials are stored in the system.	
Error message: "Invalid email format. Please enter a valid email address."	User remains on the registration page. No account is created.	
Error message: "Passwords do not match."	User remains on the registration page. No account is created.	
Mealtime and carb-to-unit ratio are updated successfully. Previous values are overwritten	User is returned to the main screen. Updated values are stored in the system.	
Error message: "Carb-to-unit ratio cannot be empty."	User remains on the edit screen. No changes are saved.	
Mealtime remains as "00:00 AM/PM" (indicating no mealtime selected). Carb-to-unit ratio is updated successfully.	User is returned to the main screen. Updated values are stored in the system.	
Mealtimes and carb-to-unit ratios for breakfast and lunch are updated successfully. Previous values are overwritten.	User is returned to the main screen. Updated values are stored in the system.	
No changes are saved. Previous values remain unchanged.	User is returned to the main screen. No updates are made to the system.	
User is authenticated, and adherence data is retrieved.	User session is initiated.	
Authentication fails. Error message is shown.	User stays on login screen.	
Error message: "Invalid username or password."	Authentication fails.	
Either session timeout on first device, or second login attempt is blocked.	One active session is maintained.	
Error message: "Username and password are required."	Authentication fails.	
Login is successful.	User session begins and adherence data is created.	
Dashboard loads with user's adherence data.	Adherence data is accessible and visible.	
Data is accessed and displayed	User is accessing data	
User is prompted to start tracking data	User begins tracking data	
User is back on the main menu	User is on the main screen	
