

Project Name:		Medicine Adherence System							
Module Name:		1.2 Authenticate User / User Log In							
Reference Document:		CMSC 355 Requirments (Use Cases)							
Created By:		Inayah O'Neil							
Date of Creation:		3/18/25							
Date of Review:									
TEST CASE ID	TEST SCENARIO	TEST CASE	PRE-CONDTION	TEST STEPS	TEST DATA	EXPECTED RESULT	POST CONDITION	ACTUAL RESULT	STATUS (PASS/FAIL)
TC_LOGIN_001	Verify credentials for user to log in to software	Enter valid User Name and valid Password	1. User has already established a valid login	1. Click "Login" Button		Successful user creation	Home page is shown		
				2. Enter User Name	<Valid Username>				
				3. Enter Password	<Valid Password>				
				4. Click "Create Login"					
TC_LOGIN_002	Verify credentials for user to log in to software	Enter invalid User Name and valid Password	1. User has already established a valid login	1. Click "Login" Button		A message showing "the username you input is invalid"			
				2. Enter User Name	<Invalid Username>				
				3. Enter Password	<Valid Password>				
				4. Click "Create Login"					
TC_LOGIN_003	Verify credentials for user to log in to software	Enter valid User Name and invalid Password	1. User has already established a valid login	1. Click "Login" Button		A message showing "the password you entered is invlaid"			
				2. Enter User Name	<Valid Username>				
				3. Enter Password	<Invalid Password>				
				4. Click "Create Login"					
TC_LOGIN_004	Verify credentials for user to log in to software	Enter invalid User Name and invalid Password	1. User has already established a valid login	1. Click "Login" Button		A message showing "the username and password you entered is invalid"			
				2. Enter User Name	<Invalid Username>				
				3. Enter Password	<Invalid Password>				
				4. Click "Create Login"					

[illegible]

[illegible]

[illegible]

Project Name:		Medicine Adherence System							
Module Name:		4.3/4.5 Send Medication Reminder							
Reference Document:		CMSC 355 Requirments (Use Cases)							
Created By:		Cara Routh							
Date of Creation:		3/25/25							
Date of Review:									
TEST CASE ID	TEST SCENARIO	TEST CASE	PRE-CONDTION	TEST STEPS	TEST DATA	EXPECTED RESULT	POST CONDITION	ACTUAL RESULT	STATUS (PASS/FAIL)
TC_4.3_4.5_001	Verify medication reminder sent and correctly handles a missed response.	Successful reminder sent and patient responds.	Patient has a medication schedule and system has contact details	1. System retrieves scheduled reminder time 2. Reminder message is generated 3. System selects the patients preffered contact method 4. The system sends the reminder 5. The systems logs the reminder as "Sent". 6. System wait for patients response within the decided timeframe. 7. System logs the response in the medicine log.	Patient ID : 12345 Contact Method : App Notification Medication : Xanax Scheduled Time : 8:00 am Response : Taken	System logs the response as taken.	The system completes tracking for this reminder.		
TC_4.3_4.5_002	Verify that the system retries sending the reminder if the patient does not respond.	Reminder is sent but patient does not respond (their dose was missed)	Patient has a medication schedule and patient doesn't respond to first reminder.	1. Send reminder. 2. The system waits for the patient's response within a set time. 3. No response is recieved 4. The system incremnets the missed count. 5. If retires are available the system resends the reminder. 6. If the retry limit is reached the system lets the caregiver know.	Patient ID : 24680 Medication : Insulin Scheduled Time : 9:00 am Retry Limit : 2	The system will resend the reminder if within the retry limit, and if the retry limit is reached, the system alerts the caregiver.	Either logs a response, retires, or escalates the situation based on patients response.		
TC_4.3_4.5_003	Verify that the system is able to handle a failure to send a reminer and retries.	Reminder failed to send because of a network issue.	Patient has a medication schedule and system attempts to send a reminder but fails.	1. Send reminder 2. Failure to send due to network issues 3. Retries sending after a set interval 4. If the retries also fail, the system logs an error	Patient ID : 56789 Medication : Percocet Failure Cause : Network Cause	The system retires sending the reminder and if all the retries fail then admin is alerted.	Either reminder is succesfully sent or an error is logged.		

[illegible]

[illegible]

[illegible]

